


Vector™

abducted shoulder sling



item #	size	measurement	
(left) 183702	small	up to 13 1/2" (34.5cm)	
183704	medium	up to 14 1/2" (37.5cm)	
183705	large	up to 15 1/2" (39.5cm)	
(right) 183712	small	up to 13 1/2" (34.5cm)	
183714	medium	up to 14 1/2" (37.5cm)	
183715	large	up to 15 1/2" (39.5cm)	

* measure from olecranon process to 5th MP joint



Medical Specialties, Inc.
4600-K Lebanon Rd.
Charlotte, NC 28227 USA
p: 800.334.4143 f: 800.694.9060
p: 704.694.2434 f: 704.694.9060
e-mail: request@medspec.com
www.medspec.com

The initial fitting of this product should be done by a medical professional. Once the brace is fit properly, the patient can remove and apply the product without having to readjust it.

Fitting Instructions for Right Shoulder


- a. Remove short strap off of abduction pillow.
- b. Unbuckle waist strap.
- c. Place abduction pillow on right side, with sling strap over left shoulder.
- d. Re-buckle waist strap and adjust strap to securely hold the pillow against body in the desired position.
- e. Place right arm in sling with elbow as far back as possible. Attach short strap at elbow to cover sling pocket.
- f. Adjust sling straps, cutting excess material if needed and then reattaching Y-tabs.
- g. Finally, adjust Vector strap to effectively reduce pressure of the sling strap on the neck.



front view with pad

Fitting Instructions for Left Shoulder


- a. Remove short strap off abduction pillow.
- b. Undo Vector strap from metal ring and flip Y-tab.
- c. Remove sling from abduction pillow leaving Velcro material on pad.
- d. Place abduction pillow at waistline of left side. Place the waist strap around the patient. Position the pad for desired internal or external rotation and adjust straps to securely hold pad in place. For excess material, strap may be cut and Y-tabs reattached.
- e. Put arm in sling, positioning the elbow as far back as possible. Attach the sling to the outside of the pillow, along the Velcro strap. Close sling and place short strap at elbow to hold in place. Bring shoulder strap over opposite shoulder. Take front shoulder strap off of D-rings. Flip Y-tab at end of shoulder strap. Re-loop through double D-rings and secure. Adjust sling straps. When properly positioned, the metal ring should be just medial to the affected side's scapula.
- f. Remove small piece of Velcro on the back of the waist strap. Rotate the Vector strap. Reattach the small piece of Velcro. Feed the Vector strap through the metal ring, pulling down on metal ring to relieve pressure on neck.
- g. Adjust all straps and pillow to desired position. The straps may be trimmed by removing the Y-tabs, cutting the strap and reattaching the Y-tab.



Diagrams:



Y-tab feature



abduction pillow

abdomen

hand

upper arm notch

Use & Care:

Hand wash in cold water using mild soap, rinse thoroughly. AIR DRY. Note: If not rinsed thoroughly, residual soap may cause irritation and damage material.