

Ryno Lacer[®]


Universal Wrist & Thumb Support



The Ryno Lacer provides superior immobilization of the CMC and MP joints of the thumb, yet it permits limited motion of the IP joint to allow for contact with the opposing fingers. The universal design fits a broad range of patient sizes and minimizes inventory costs.

Indications:

- DeQuervain's Tendonitis
- Gamekeeper's Thumb
- Post Cast Healing
- Selected Scaphoid Injuries



Universal Design

Dorsal stay casing may easily be repositioned to fit a wide range of patient sizes.

3 Stay Stabilization

Radial, palmar, and dorsal aluminum stays immobilize the wrist and thumb. Palmar stay may be removed if desired.

Ease of Application

Lace closure allows for easy one-handed application.

Optimum Comfort

The polypropylene felt liner and outer nylon knit fabric allow moisture to evaporate away from the skin and enhance patient compliance.

Suggested HCPCS: L3807

long left	long right	size	wrist circ.
223520	223521	univ.	5 1/2" - 8 1/2"

short left	short right	size	wrist circ.
223510	223511	univ.	5 1/2" - 8 1/2"

