

# CMC-X™ Lacer

## Thumb Stabilizer



The CMC-X™ Lacer Thumb Stabilizer provides your patients with pain relief from carpometacarpal (CMC) joint arthritis while permitting them to do the activities they enjoy.

### “X” Strapping System

The X straps are easily adjusted to position the first metacarpal (base of thumb) in the desired degree of abduction to help stabilize the CMC joint while still allowing the thumb to have functional mobility.

### Adjustable Stay Pod

The adjustable stay pod works in concert with the “X” straps to provide gentle pressure and help stabilize the CMC joint. The stay pod contains a malleable aluminum stay and is easily bent and repositioned by a medical provider to provide a customized fit.

### Lacer Closure System

The lacer closure system enhances ease of application.

### Bilateral Design

The bilateral design fits the left or right thumb which reduces the carrying cost of inventory.

### Indications:

- CMC joint osteoarthritis (OA)
- Medial subluxation of the CMC joint



Size	Item #	Measurements
x-small	224471	under 5½”
small	224472	5½” - 6¼”
medium	224474	6¼” - 7”
large	224475	7” - 7¾”
x-large	224476	over 7¾”

\*measure circumference of wrist at narrowest point

Suggested HCPCS: L3923 or L3924