

EVO® Quatro® Ankle Stabilizer



The next EVolution in ankle bracing technology™

The EVO Quatro provides the innovative technology of the EVO's new strapping and cuff system in a low-profile design, making it the ideal stabilizer for those athletes or patients that will not tolerate any bulky materials inside their shoe, yet need strong ankle support.



Low-profile Design

The boot design of the Quatro has been modified to reduce the coverage around the mid-foot and allow for unrestricted plantar flexion.

Stirrup & Stabilizing Straps

Internal stirrup strap travels through slots in the dynamic cuff and continues outside the boot to form the stabilizing straps. This patented design better captures the calcaneus (heel) and resists inversion or eversion of the ankle to help stabilize the talocrural and subtalar joints.

Dynamic Cuff

Ultra thin plastic cuff works in concert with strapping system to provide circumferential compression around the distal tibiofibular joint. In addition, the dynamic cuff is now lined with CoolFlex™ padding for improved comfort and Skinloc™ material to resist distal migration.



Elastic cuff has been removed for illustrative purposes



Indications:

- Inversion or eversion ankle sprains
- Syndesmotic ankle sprains (High ankle sprains)
- Medial or lateral ankle instability
- Subtalar joint instability

Bilateral design fits left or right foot

Black	White	Size	Ankle Circum.*
264231	264221	x-small	10" - 11"
264232	264222	small	11" - 12"
264234	264224	medium	12" - 13"
264235	264225	large	13" - 14"
264236	264226	x-large	14" - 15"
264237	264227	xx-large	15" - 16"

*measure ankle circumference according to diagram below



Suggested HCPCS: L1902

Men's Shoe Size	6-7	7 1/2-9	9 1/2-11	11 1/2-13	13 1/2-15
Narrow Foot		xs	s	m	l
Regular Foot	xs	s	m	l	xl
Wide Foot	s	m	l	xl	xxl
Women's Shoe Size	5-6	6 1/2-8	8 1/2-10	10 1/2-12	12 1/2-14
Narrow Foot		xs	s	m	l
Regular Foot	xs	s	m	l	xl
Wide Foot	s	m	l	l	xl

Medical Specialties, Inc.

4600-K Lebanon Rd. Charlotte, NC 28227 p: 800-334-4143 f: 800-694-9060 p: 704-694-2434 f: 704-694-9060

e-mail: request@medspec.com www.medspec.com

Patent Numbers: US 8,721,578; US 8,808,215; US 9,375,339; US 9,737,430; US 10,383,404; AU 2008317265; CA 2,703,266;

European (AT, BE, CH, DE, DK, FR, GB, IE, NL, NO, SE) 2203135; patent pending

ASO and EVO are registered trademarks of Medical Specialties, Inc.

© 2014-2021 Medical Specialties, Inc.

