

Supinator™

PTT Stabilizer



The Supinator is a comfortable and relatively low-profile brace for the treatment of posterior tibial tendon dysfunction (PTTD). With its patent pending design, the Supinator reduces over-pronation by applying a supinating force.

Arch Strap & Pad

The arch strap works in concert with a repositionable arch pad to lift the medial arch. The arch pad may also be removed if desired and to allow the Supinator to be worn with an orthotic.

Heel Strap

The heel strap restricts eversion of the calcaneus (heel).

Speed Lacer

Speed Lacer II closure system allows for ease of application and provides equal tension across all laces.

Contraindications:


- Lateral ankle sprain
- Lateral ankle instability

Indications:

- Posterior tibial tendon dysfunction (PTTD)
- Acquired adult flatfoot deformity
- Posterior tibial tendon insufficiency



left foot pictured above



right foot pictured above



left foot pictured above

Suggested HCPCS: L1902

Left	Right	Size	Ankle Circum.*
264251	264261	X-Small	10" - 11"
264252	264262	Small	11" - 12"
264254	264264	Medium	12" - 13"
264255	264265	Large	13" - 14"
264256	264266	X-Large	14" - 15"

*measure ankle circumference according to diagram below



Men's Shoe Size	6-7	7 1/2-9	9 1/2-11	11 1/2-13	13 1/2-15
Narrow Foot		xs	s	m	l
Regular Foot	xs	s	m	l	xl
Wide Foot	s	m	l	xl	
Women's Shoe Size	5-6	6 1/2-8	8 1/2-10	10 1/2-12	12 1/2-14
Narrow Foot		xs	s	m	l
Regular Foot	xs	s	m	l	xl
Wide Foot	s	m	l	xl	

Medical Specialties, Inc. 4600-K Lebanon Rd. Charlotte, NC 28227
 p: 800-334-4143 f: 800-694-9060 p: 704-694-2434 f: 704-694-9060
 e-mail: request@medspec.com www.medspec.com Patent Pending
 Supinator is a trademark of Medical Specialties, Inc. © 2018-2021 Medical Specialties, Inc.



rev 9/21