

Archimed® 63I

Universal Spinal Brace



The Archimed 63I Universal Spinal Brace's patented design easily adjusts to fit patients of various sizes and provides the support they need in the treatment of lumbar spine disorders.



Universal Design

Attached outer panels easily slide to fit patients with waist sizes ranging from 29" to 57". An available extension panel fits patients with a waist size up to 72". Both panels are now more breathable and lower profile for enhanced comfort.



Rigid Anterior and Posterior Panels

Polyethylene plastic panels work in concert with the pulley closure system to easily generate radial compression and increase intra-abdominal pressure. The posterior panel extends from sacrococcygeal Junction to T-9 vertebra.



Indications:

- Post-operative lumbar spine stabilization
- Degenerative disc disease
- Ligamentous sprains & muscle strains of the lumbar spine

PDAC Verified HCPCS:
L063I or L0648

New Pulley Closure System

New cord dual pulley closure system is easier to pull and provides a combined 6:1 mechanical advantage. The durable cord pulley closure system will never come off track.

Product	Size	Item #	Waist Circ.	PDAC Verified
Archimed 63I	Universal	163310	29" - 57"	L063I or L0648
Archimed Extension Panel		163390	57" - 72"	



Medical Specialties, Inc. 4600-K Lebanon Rd. Charlotte, NC 28227 p: 800-334-4143 f: 800-694-9060
p: 704-694-2434 f: 704-694-9060 e-mail: request@medspec.com www.medspec.com US Patent # 9,808,369
©2014-2020 Medical Specialties, Inc.

rev 12/20