

Archimed[®] 627

Universal Spinal Brace



The Archimed 627 Universal Spinal Brace's patented design easily adjusts to fit patients of various sizes and provides the support they need in the treatment of lumbar spine disorders.

Universal Design

Attached outer panels easily slide to fit patients with waist sizes ranging from 29" to 57". An available extension panel fits patients with a waist size up to 72". Both panels are now more breathable and lower profile for enhanced comfort.



Rigid Anterior and Posterior Panels

Polyethylene plastic panels work in concert with the pulley closure system to easily generate radial compression and increase intra-abdominal pressure. The posterior panel extends from L-1 to below L-5 vertebra.



Indications:

- Post-operative lumbar spine stabilization
- Degenerative disc disease
- Ligamentous sprains & muscle strains of the lumbar spine

**PDAC Verified HCPCS:
L0627 or L0642**

New Pulley Closure System

New cord dual pulley closure system is easier to pull and provides a combined 6:1 mechanical advantage. The durable cord pulley closure system will never come off track.

Product	Size	Item #	Waist Circ.	PDAC Verified
Archimed 627	Universal	163320	29" - 57"	L0627 or L0642
Archimed Extension Panel		163390	57" - 72"	

