

# Gripper™

## 12" Hinged Knee Brace



The Gripper™ Hinged Knee Brace provides your athletes and patients with the medial and lateral stability they need while resisting distal migration.



### Compliance & Migration

Internal SkinLoc™ pad works in concert with the outer closure strap to grip the proximal end of the calf muscle and inhibit the brace from migrating distally, thereby enhancing patient compliance

### Durability

Straps are made of high tension elastic that allows for numerous applications without fraying or losing closure strength. Metal hinge pockets are reinforced on the inside with ballistic nylon material at wear points to prevent hinge from wearing through the pocket.

### Ease of Application & Comfort

The wraparound anterior closure allows for ease of application and an optimum fit. Open popliteal minimizes material bunching during flexion to improve comfort.



The Gripper 12" is available in two different base materials: Neoprene and CoolFlex™. The neoprene version provides warmth and superior compression. The Gripper 2 is made with a thinner and lighter neoprene. The CoolFlex version is breathable and hypoallergenic.

### Indications:

- Meniscal injuries
- MCL or LCL injuries or sprains
- Medial or lateral knee instabilities

Suggested HCPCS: L1820

Gripper Neoprene 3/16"	Gripper CoolFlex	Gripper 2 Neoprene 1/8"	Size	6" Above	Mid-Patella	6" Below
117011	117031	117151	X-Small	11" - 13"	10" - 12"	10" - 12"
117012	117032	117152	Small	13" - 16"	12" - 14"	12" - 14"
117014	117034	117154	Medium	16" - 18"	14" - 16"	13" - 15"
117015	117035	117155	Large	18" - 21"	16" - 18"	15" - 17"
117016	117036	117156	X-Large	21" - 24"	18" - 20"	17" - 19"
117017	117037	117157	XX-Large	24" - 27"	20" - 22"	19" - 21"
117018	117038	117158	3X-Large	27" - 31"	22" - 24"	21" - 23"
117019	117039	117159	4X-Large	31" - 35"	24" - 26"	23" - 26"
117010	117030		5X-Large	35" - 39"	26" - 28"	26" - 29"

\*Use mid-patella measurement for primary sizing info

