

Gripper™

16" Hinged Knee Brace



The Gripper™ Hinged Knee Brace provides your athletes and patients with the medial and lateral stability they need while resisting distal migration.



Compliance & Migration

Internal SkinLoc™ pad works in concert with the outer closure strap to grip the proximal end of the calf muscle and inhibit the brace from migrating distally, thereby enhancing patient compliance

Durability

Straps are made of high tension elastic that allows for numerous applications without fraying or losing closure strength. Metal hinge pockets are reinforced on the inside with ballistic nylon material at wear points to prevent hinge from wearing through the pocket.

Ease of Application & Comfort

The wraparound anterior closure allows for ease of application and an optimum fit. Open popliteal minimizes material bunching during flexion to improve comfort.

Indications:

- Meniscal injuries
- MCL or LCL injuries or sprains
- Medial or lateral knee instabilities

Suggested HCPCS: L1820



geared polycentric

aluminum

The Gripper 16" is available in two different base materials: Neoprene and CoolFlex™. The neoprene version provides warmth and superior compression. The CoolFlex version is breathable and hypoallergenic.

Neoprene Geared	Neoprene Aluminum	CoolFlex Geared	CoolFlex Aluminum	Size	6" Above	Mid-Patella	6" Below
117072	117072-A	117082	117082-A	Small	13" - 16"	12" - 14"	12" - 14"
117074	117074-A	117084	117084-A	Medium	16" - 18"	14" - 16"	13" - 15"
117075	117075-A	117085	117085-A	Large	18" - 21"	16" - 18"	15" - 17"
117076	117076-A	117086	117086-A	X-Large	21" - 24"	18" - 20"	17" - 19"
117077	117077-A	117087	117087-A	XX-Large	24" - 27"	20" - 22"	19" - 21"

*Use mid-patella measurement for primary sizing info



rev 4/23

Medical Specialties, Inc. 4600-K Lebanon Rd. Charlotte, NC 28227 p: 800-334-4143 f: 800-694-9060
 p: 704-694-2434 f: 704-694-9060 e-mail: request@medspec.com www.medspec.com
 ©2010-2023 Medical Specialties, Inc.

