

Ryno Lacer® II

Universal Wrist & Thumb Support



With its patent pending design, the new Ryno Lacer II Universal provides equal tension across all laces offering supreme comfort, ease of application and superior immobilization of the wrist and thumb.

Lacer II Closure System

Unique lace closure design provides equal tension across all laces to easily generate uniform compression for unparalleled support. The tension can be easily adjusted for patients with disproportionate forearm and wrist circumferences.

Universal Design

Dorsal stay casing can be easily repositioned to fit a wide range of patient sizes. In addition, dorsal, palmar and radial stays are removable if desired.

Comfort and Compliance

New softer proprietary felt liner contains polypropylene and wicks moisture away from the skin. Memory foam padding in palm and around base of thumb enhances comfort and support.

Indications:

- DeQuervain's Tendonitis
- Gamekeeper's Thumb
- Post Cast Healing
- Select Scaphoid Injuries

Long Left	Long Right	Size	Wrist Circ.
223760	223761	univ.	5 1/2" - 8 1/2"
Short Left	Short Right	Size	Wrist Circ.
223770	223771	univ.	5 1/2" - 8 1/2"

Suggested HCPCS: L3807 or L3809



Long Ryno Lacer II



Short Ryno Lacer II



Short Ryno Lacer II



rev 6/23

Medical Specialties, Inc. 4600-K Lebanon Rd. Charlotte, NC 28227 p: 800-334-4143 f: 800-694-9060 p: 704-694-2434 f: 704-694-9060

e-mail: request@medspec.com www.medspec.com US Patent No. 11,660,221

Ryno Lacer is a registered trademark of Medical Specialties, Inc. © 2017-2023 Medical Specialties, Inc.

