

Archimed® 637

Spinal Brace



The Archimed 637 Spinal Brace utilizes a compound pulley system to comfortably and easily provide your patients with the support and stability they need in the treatment of injuries and diseases of the lower spine.

Rigid Anterior, Lateral and Posterior Panels

The anterior, lateral and posterior panels work in concert with the pulley-lace closure system to easily generate intra-abdominal pressure and restrict coronal and sagittal plane motion. The posterior panel extends from the sacrococcygeal junction to the T-9 vertebra with a central cut-out over the lumbar region for post operative application. Posterior panel and interior padding also contour automatically to the patient's lordotic curve to provide comfortable support.

Pulley-Lace Closure System

The dual pulley-lace closure system provides a combined 6 to 1 mechanical advantage. In addition, the pulley-lace closure system is adjustable to provide a wide range of fit per given size and automatically conforms to various body shapes. The durable pulley-lace closure system will never come off track.



PDAC Verified HCPCS: L0637 or L0650

Indications:

- Post-operative lumbar spine stabilization
- Degenerative disc disease
- Ligamentous sprains & muscle strains of the lumbar spine



Size	Item #	Waist Circ.
Small	163332	27" - 31"
Medium	163334	32" - 38"
Large	163335	39" - 45"
X-Large	163336	46" - 52"
XX-Large	163337	53" - 59"
3X-Large	163338	over 59"