

ASO[®] with Plastic Stays Ankle Stabilizer



The Choice of ChampionsSM

With the new improvements to the ASO with Stays, the ASO maintains its leadership position in the ankle brace market. Give your athletes and patients the patented and proven technology of the ASO with Stays. The leader continues to stay several steps ahead...

Plastic Stays

Removable plastic stays allow for enhanced proprioceptive support.

Stabilizing Straps

The stabilizing straps are now angled to better contour above the malleoli for enhanced comfort and better migration resistance. Nylon stabilizing straps form a complete figure-eight to help stabilize the subtalar and talocrural joints.

Elastic Cuff

Circumferential woven elastic cuff keeps laces and stabilizing straps secure. In addition, it works in concert with the upper boot to resist the inversion or eversion forces.

CoolFlex Achilles Pad & Tongue

CoolFlexTM padding has been increased on the inside of the boot along the Achilles tendon to provide greater comfort. The shoe lace is attached to the tongue to ensure lace is always centered.

Invisible Seam

There is no longer a seam under the foot. The seam is now hidden along the rear portion of the boot to improve comfort.



Indications:

- Inversion or Eversion Ankle Sprains
- Medial or Lateral Ankle Instability
- Subtalar Joint Instability

PDAC Verified HCPCS:
L1902

| Men's Shoe Size | 6-7 | 7 1/2-9 | 9 1/2-11 | 11 1/2-13 | 13 1/2-15 |
|-------------------|-----|---------|----------|-----------|-----------|
| Narrow Foot | | xs | s | m | l |
| Regular Foot | xs | s | m | l | xl |
| Wide Foot | s | m | l | xl | xxl |
| Women's Shoe Size | 5-6 | 6 1/2-8 | 8 1/2-10 | 10 1/2-12 | 12 1/2-14 |
| Narrow Foot | | xs | s | m | l |
| Regular Foot | xs | s | m | m | l |
| Wide Foot | s | m | l | l | xl |

Bilateral design fits left or right foot

| Black | White | Size | Ankle Circum.* |
|--------|--------|----------|----------------|
| 264031 | 223181 | x-small | 10" - 11" |
| 264032 | 223182 | small | 11" - 12" |
| 264034 | 223184 | medium | 12" - 13" |
| 264035 | 223185 | large | 13" - 14" |
| 264036 | 223186 | x-large | 14" - 15" |
| 264037 | 223187 | xx-large | 15" - 16" |

*measure ankle circumference according to diagram below



Medical Specialties Inc. 4600-K Lebanon Rd. Charlotte, NC 28227
p: 800-334-4143 f: 800-694-9060 p: 704-694-2434 f: 704-694-9060

e-mail: request@medspec.com www.medspec.com

Patents: US 8,808,215; US 9,375,339; US 10,383,404 and other patents pending

ASO is a registered trademark of Medical Specialties Inc. © 2009-2024 Medical Specialties, Inc.