

EVO[®] Hinge Ankle Stabilizer



The EVO Hinge Ankle Stabilizer is an evolutionary step forward in ankle protection. With its patented design and removable hinge, the EVO Hinge provides even greater inversion/eversion control and protection against syndesmotic ankle sprains than the ASO Flex-Hinge.



actual product

Removable Plastic Hinge

New lower profile plastic hinge fits in almost any shoe style. Removable plastic hinge allows for a progressive step down treatment protocol.

Dynamic Cuff

Internal plastic cuff applies circumferential compression around the distal tibiofibular joint. The dynamic cuff increases the circumferential compression as the stabilizing straps restrict inversion or eversion of the calcaneus (heel).

Stirrup Strap

Inner stirrup strap is lined with Skinloc[™] material and connected to internal plastic cuff to better capture the calcaneus (heel) and resist inversion or eversion of the ankle.

Stabilizing Straps

Nylon stabilizing straps form a figure-eight to help stabilize the talocrural and subtalar joints.

Indications:

- Inversion or eversion ankle sprains
- Syndesmotic ankle sprains (high ankle sprains)
- Medial or lateral ankle instability
- Subtalar joint instability
- Midtarsal joint instability



The next evolution in ankle bracing technology[™]

Bilateral design fits left or right foot

Black	Size	Ankle Circum.*
264311	x-small	10" - 11"
264312	small	11" - 12"
264314	medium	12" - 13"
264315	large	13" - 14"
264316	x-large	14" - 15"
264317	xx-large	15" - 16"

PDAC Verified HCPCS:
L1906

Men's Shoe Size	6-7	7 1/2-9	9 1/2-11	11 1/2-13	13 1/2-15
Narrow Foot		xs	s	m	l
Regular Foot	xs	s	m	l	xl
Wide Foot	s	m	l	xl	xxl
Women's Shoe Size	5-6	6 1/2-8	8 1/2-10	10 1/2-12	12 1/2-14
Narrow Foot		xs	s	m	l
Regular Foot	xs	s	m	l	xl
Wide Foot	s	m	l	l	xl

*measure ankle circumference according to diagram below



Medical Specialties Inc.
4600-K Lebanon Rd. Charlotte, NC 28227 p: 800-334-4143 f: 800-694-9060 p: 704-694-2434 f: 704-694-9060
e-mail: request@medspec.com www.medspec.com
Patent Numbers: US 8,721,578; US 8,808,215; US 9,375,339; US 9,737,430; US 10,383,404; AU 2008317265; CA 2,703,266;
European (AT, BE, CH, DE, DK, FR, GB, IE, NL, NO, SE) 2203135; patent pending
ASO and EVO are registered trademarks of Medical Specialties Inc.
© 2013-2024 Medical Specialties, Inc.

