

AKS™ Knee Support with Plastic Hinges



The AKS Knee Support offers your patients the unique ability to adjust compression and fit from top to bottom.



Adjustable Knee Sleeve

The distal and proximal openings allow for easy application and ability to fit various patient shapes. In addition, compression can be adjusted to the desired level.

Adjustable Patellar Buttress

Soft horseshoe buttress is repositionable and works in concert with central sleeve region to provide circumferential compression around the knee and focal compression around the patella.

Medial and Lateral Plastic Hinges

Plastic hinges enhance the brace's ability to conform to the knee during flexion. Internal Skinloc™ pad and plastic hinges function synergistically to resist distal migration and vertical bunching.



Indications:

- Arthritis of the knee
- Bursitis of the knee
- Chondromalacia patella
- Patellar tendinitis

Bilateral design fits left or right knee.
Neoprene version provides greater compression and warmth.
CoolFlex version is breathable and lightweight for superior comfort.

Neoprene	Size	Knee Circ.*	CoolFlex
117091	X-Small	10" - 12"	117101
117092	Small	12" - 14"	117102
117094	Medium	14" - 16"	117104
117095	Large	16" - 18"	117105
117096	X-Large	18" - 20"	117106
117097	XX-Large	20" - 22"	117107
117098	3X-Large	over 22"	117108

Suggested HCPCS: L1810



*measure circumference of knee at mid patella

Medical Specialties Inc.
4600-K Lebanon Rd.
Charlotte, NC 28227
p: 800-334-4143 f: 800-694-9060
p: 704-694-2434 f: 704-694-9060
e-mail: request@medspec.com
www.medspec.com
©2010-2024 Medical Specialties Inc.