

# PFO<sup>®</sup> Plantar Fascia Orthosis



**PFO<sup>®</sup> Plantar Fascia Orthosis** is a comfortable day or night treatment option for plantar fasciitis.

### Elastic Tensioner

High strength elastic strap applies uniform pressure under plantar fascia to provide continuous tension and stretching. Tension is adjustable and strap is repositionable.

### Breathable Ankle Wrap

Polypropylene lining wicks moisture away from the skin and will not support bacterial growth unlike most foam materials.

### Low Profile

Low-profile design allows PFO to easily fit inside shoe. PFO may also be worn directly next to the skin.

### Indications:

- Plantar Fasciitis
- Plantar Flexion Contracture



Bilateral design fits left or right foot.

Suggested HCPCS: L1902

Item #	Size	Womens's US Shoe Size	Men's US Shoe Size
264144	Small/Medium	up to 9.5	up to 9
264146	Large/X-Large	10+	9.5+



Medical Specialties Inc.  
4600-K Lebanon Rd. Charlotte, NC 28227 p: 800-334-4143 f: 800-694-9060 p: 704-694-2434 f: 704-694-9060  
e-mail: request@medspec.com www.medspec.com  
© 2020-2024 Medical Specialties Inc.

