

Thumb Supports



TeePee™ Thumb Protector

This unique configuration with its three aluminum stays forms the “TeePee” that comfortably supports the MCP (metacarpal-phalangeal) joint. Each stay can be shaped and contoured with ease to provide a custom fit. The TeePee does not restrict wrist motion and fits the left or right thumb (bilateral).



PDAC Verified HCPCS: L3923

Sizing for the TeePee
To Determine proper size, measure the linear distance from the dorsal crease of the thumb IP joint to the peak of the radial styloid process.



| Size | Item Number | Measurements |
|---------|-------------|-----------------|
| x-small | 223081 | under 3 1/2" |
| small | 223082 | 3 1/2" - 4 1/8" |
| medium | 223084 | 4 1/8" - 4 5/8" |
| large | 223085 | 4 5/8" - 5 1/8" |
| x-large | 223086 | over 5 1/8" |

CMC Thumb Support

The CMC Thumb Support aids in the treatment of arthritis and instability of the CMC (carpometacarpal) joint. The unique design of the internal padded stay in combination with the outer elastic strap provides focal compression at the base of the thumb and helps stabilize the CMC joint. The CMC Thumb Support provides this functional stabilization without restricting the use of the thumb in pinching or grasping objects.



Suggested HCPCS: L3923

| Size | Left | Right | Measurements |
|---------|--------|--------|-----------------|
| x-small | 224421 | 224431 | under 5 1/2" |
| small | 224422 | 224432 | 5 1/2" - 6 1/4" |
| medium | 224424 | 224434 | 6 1/4" - 7" |
| large | 224425 | 224435 | 7" - 7 3/4" |
| x-large | 224426 | 224436 | over 7 3/4" |

*measure circumference of wrist at narrowest point

Medical Specialties Inc.
4600-K Lebanon Rd.
Charlotte, NC 28227
p: 800-334-4143 f: 800-694-9060
p: 704-694-2434 f: 704-694-9060
e-mail: request@medspec.com
www.medspec.com



©2007-2024 Medical Specialties Inc.



rev 3/24