

# Viper™ Wrist Support



The Viper provides comfortable compression and support while moderately restricting wrist motion. The ease of application and low profile design ensures a high rate of patient compliance.

## Indications:

- Carpal tunnel syndrome
- Wrist sprains or strains
- Post cast removal
- Post operative care

## Optimum Comfort

The Viper Wrist Support is made from our unique CoolFlex™ material that is soft, breathable, and provides uniform compression.

## Restricted Motion

The aluminum palmar stay resists flexion of the wrist and is malleable to allow for a custom fit. There is no dorsal stay to resist extension of the wrist.

## V-Strap Retention System

Allows extraordinary ease of application since the “V” formed by the compression strap and the hook fastener prevents the strap from sliding out of the internal slot when applying or removing the wrist support.

Suggested HCPCS: L3908

Universal Left	Universal Right	Wrist Circ.
223530	223531	4 1/2" - 9"

\*measure circumference of wrist at narrowest point



Medical Specialties Inc. 4600-K Lebanon Rd. Charlotte, NC 28227 p: 800-334-4143 f: 800-694-9060 p: 704-694-2434 f: 704-694-9060  
e-mail: request@medspec.com www.medspec.com  
©2008-2024 Medical Specialties Inc.

