



# Wrist Lacer™

## wrist support

With the ease of application of the Med Spec Wrist Lacer, you don't have to worry about your patients being compliant. In addition, they will love the exceptional comfort and support the Wrist Lacer provides.



8" wrist support

### Ease of Application

The lace closure mechanism provides unparalleled ease of application. You only need one free hand to tension and fasten the lace closure. The stockinette sleeve prevents the thumb from catching the laces and provides padding between the laces and the patient's skin.

### Adjustability & Immobilization

Malleable palmar and dorsal aluminum stays may be bent to the desired degree of flexion or extension. The blue palmar pad provides enhanced support by matching the contour of the depression in the palm.

### Optimum Comfort

The anatomically contoured design and the adjustable thumb strap provide a comfortable fit and allow unrestricted motion of the fingers.

### Indications:

- Carpal tunnel syndrome
- Wrist sprains or strains
- Post cast removal
- Post operative care



10 1/2" wrist & forearm support

Now available in Pediatric (xx-small) size.

Suggested HCPCS: L3908

Size	8" black left	8" black right
xx-small	223929	223939
x-small	223921	223931
small	223922	223932
medium	223924	223934
large	223925	223935
x-large	223926	223936

Sizing Chart	
Size	Wrist Circ.
xx-small	under 4 3/4"
x-small	4 3/4" - 5 1/2"
small	5 1/2" - 6 1/4"
medium	6 1/4" - 7"
large	7" - 7 3/4"
x-large	7 3/4" - 8 1/2"

Size	10 1/2" black left	10 1/2" black right
small	223952	223962
medium	223954	223964
large	223955	223965
x-large	223956	223966

