

Wrist Lacer™ II

wrist support



With its patented design, the Wrist Lacer II provides equal tension across all laces offering supreme comfort, immobilization and ease of application.

Adjustability & Immobilization

Blue palmar pad offers enhanced support by matching the contour of the palm. The dorsal and palmar aluminum stays are removable to allow bending by medical professionals to create a custom fit.

Lacer II Closure System

Unique lace closure design provides equal tension across all laces to easily generate uniform compression for unparalleled support. The tension is easily adjusted for patients with unusual forearm and wrist circumferences.

Optimum Comfort

Pad inside stockinette protects tendons by cushioning distal lace at base of thumb. The polypropylene felt liner and perforated suede allow moisture to evaporate away from the skin thereby enhancing patient compliance.



8" wrist support



10 1/2" wrist & forearm support

Indications:

- Carpal tunnel syndrome
- Wrist sprains or strains
- Post cast removal
- Postoperative care

Suggested HCPCS: L3908

| Size | 8" Back Left | 8" Black Right |
|-----------|--------------|----------------|
| Pediatric | 223329 | 223339 |
| X-Small | 223321 | 223331 |
| Small | 223322 | 223332 |
| Medium | 223324 | 223334 |
| Large | 223325 | 223335 |
| X-Large | 223326 | 223336 |
| XX-Large | 223327 | 223337 |

| Sizing Chart | |
|--------------|-----------------|
| Size | Wrist Circ. |
| xx-small | under 4 3/4" |
| x-small | 4 3/4" - 5 1/2" |
| small | 5 1/2" - 6 1/4" |
| medium | 6 1/4" - 7" |
| large | 7" - 7 3/4" |
| x-large | 7 3/4" - 8 1/2" |
| xx-large | 8 1/2" - 10" |

| Size | 10 1/2" Black Left | 10 1/2" Black Right |
|---------|--------------------|---------------------|
| x-small | 223361 | 223371 |
| small | 223362 | 223372 |
| medium | 223364 | 223374 |
| large | 223365 | 223375 |
| x-large | 223366 | 223376 |

