

Wrist Lacer™

Universal Wrist Support



With the ease of application and exceptional comfort of the Med Spec Wrist Lacer, you don't have to worry about your patients being compliant. In addition, the universal design fits a broad range of patient sizes and minimizes inventory costs.



Indications:

- Carpal tunnel syndrome
- Wrist sprains or strains
- Post cast removal
- Post operative care

Universal Design

Dorsal stay casing may easily be repositioned to fit a wide range of patient sizes.

Ease of Application

Lace closure allows for easy one-handed application.

Immobilization

Malleable palmar and dorsal aluminum stays may be bent to the desired degree of flexion or extension. Palmar stay may be removed if desired.

Optimum Comfort

The polymer felt liner and outer nylon knit fabric allow moisture to evaporate away from the skin and enhance patient compliance. The anatomically contoured design and the adjustable thumb strap provide a comfortable fit and allow unrestricted motion of the fingers and thumb.



8" left	8" right	size	wrist circ.
223490	223491	univ.	5 ½" - 8 ½"
10 ½" left	10 ½" right	size	wrist circ.
223500	223501	univ.	5 ½" - 8 ½"

Suggested HCPCS: L3908



Medical Specialties Inc. 4600-K Lebanon Rd. Charlotte, NC 28227
p: 800-334-4143 f: 800-694-9060 p: 704-694-2434 f: 704-694-9060
e-mail: request@medspec.com www.medspec.com
© 2008-2024 Medical Specialties Inc.

