

Archimed[®] 627

Spinal Brace



The Archimed 627 Spinal Brace utilizes a compound pulley system to comfortably and easily provide your patients with the support and stability they need in the treatment of injuries and diseases of the lower spine.

Pulley-Lace Closure System

The dual pulley-lace closure system provides a combined 6 to 1 mechanical advantage. In addition, the pulley-lace closure system is adjustable to provide a wide range of fit per given size and automatically conforms to various body shapes. The durable pulley-lace closure system will never come off track.

Rigid Anterior and Posterior Panels

The polyethylene plastic panels work in concert with the pulley-lace closure system to easily generate radial compression and increase intra-abdominal pressure. The posterior panel extends from L-1 to above L-5 vertebra.

Indications:

- Post-operative lumbar spine stabilization
- Degenerative disc disease
- Ligamentous sprains & muscle strains of the lumbar spine



PDAC Verified HCPCS: L0627 or L0642

Size	Item #	Waist Circ.
Small	163322	27" - 31"
Medium	163324	32" - 38"
Large	163325	39" - 45"
X-Large	163326	46" - 52"
XX-Large	163327	53" - 59"
3X-Large	163328	over 59"



Medical Specialties Inc. 4600-K Lebanon Rd. Charlotte, NC 28227 p: 800-334-4143 f: 800-694-9060
 p: 704-694-2434 f: 704-694-9060 e-mail: request@medspec.com www.medspec.com

©2012-2024 Medical Specialties Inc.

