

# Back-n-Black™

## Back Support



The patented design features of the Back-n-Black™ provide your patients with the support and comfort they need in the treatment of lumbar spine related disorders.

### Anterior and Posterior Panels

The semi-rigid plastic anterior panel and plastic & foam posterior panels work in concert with the inelastic nylon strapping system to increase the intra-abdominal pressure, thereby decreasing the load on the intervertebral discs and surrounding musculoskeletal structures of the lumbar spine.

### Nylon Strapping System

The new larger handles on the elastic belt and nylon strapping system allow for ease of application and tightening to create the desired level of compression.

### Breathable Elastic Panels

The breathable elastic panels on the interior offer adjustability in fitting various abdominal shapes and enhance patient compliance with their low profile design.

### Indications:

- Chronic or acute low back pain
- Ligamentous sprains & muscle strains of the lumbar spine
- Postural support



| Item # | Size      | Waist Circumference |
|--------|-----------|---------------------|
| 163201 | x-small   | 22" - 26"           |
| 163202 | small     | 27" - 31"           |
| 163204 | medium    | 32" - 38"           |
| 163205 | large     | 39" - 45"           |
| 163206 | x-large   | 46" - 52"           |
| 163207 | xx-large  | 53" - 59"           |
| 163208 | xxx-large | over 59"            |

**PDAC Verified HCPCS: L0626 or L0641**

