

# Med Spec<sup>®</sup> Boxer Splint



For your patients with a boxer's fracture, the Med Spec Boxer Splint easily adjusts to the desired position and offers superior comfort and support to enhance compliance.

### Adjustable Design

Malleable palmar and dorsal aluminum stays are easily bent to immobilize the patient's metacarpophalangeal (MP) joints in the desired degree of flexion and the interphalangeal (IP) joints fully extended.

### Comfort & Support

The internal soft plastic clamshell protects and supports the fourth and fifth metacarpals without restricting the flexion of the MP joints. The polypropylene felt liner provides superb padding and allows moisture to rapidly evaporate away from the skin.

### Ease of Application

The metacarpal closure strap features our patented V-Strap Retention System that prevents the strap from sliding out of the plastic buckle when applying or removing the splint. In addition, the finger and wrist closure straps are made with a soft elastic material that may easily be adjusted as swelling increases or subsides.



### Indications:

- Fourth metacarpal fracture
- Fifth metacarpal fracture

Suggested HCPCS: L3807

Size	Left	Right	Wrist Circumference
X-Small	223571	223581	4" - 5"
Small	223572	223582	5" - 6"
Medium	223574	223584	6" - 7"
Large	223575	223585	7" - 8 <sup>1</sup> / <sub>4</sub> "
X-Large	223576	223586	8 <sup>1</sup> / <sub>4</sub> " - 9 <sup>1</sup> / <sub>2</sub> "



Medical Specialties Inc. 4600-K Lebanon Rd. Charlotte, NC 28227 p: 800-334-4143 f: 800-694-9060  
 p: 704-694-2434 f: 704-694-9060 e-mail: request@medspec.com www.medspec.com  
 U.S. Patent # 8,246,560 ©2009-2024 Medical Specialties Inc.

