

Compressor™

wrist support



The Compressor™ wrist support provides excellent comfort and compression along the carpal region. The patented strap retention system on the end of the Compressor strap allows for ease of application. In addition, the Compressor strap works in conjunction with the plastic stays to create a functional wrist support that limits the range of flexion and extension of the wrist.

Compressor Strap

The patented strap retention system allows for ease of application. High quality elastic strap also provides superior compression along the carpal region.

CoolFlex™ Material

The Compressor body is made from this breathable material that offers cushioning and four-way stretch.

Molded Plastic Palmar Stay

This high density polyethylene plastic stay is molded to provide an anatomical fit and resists flexion of wrist.

Plastic Dorsal Stay

This low density plastic stay resists extension of wrist.

Indications:

- Carpal tunnel syndrome
- Wrist sprains or strains
- Arthritis of the wrist

Suggested HCPCS: L3908

Left	Right	Size	Wrist Circ.
223802	223812	small	5" - 6"
223804	223814	medium	6" - 7"
223805	223815	large	7" - 8"
223806	223816	x-large	8" - 9"

