

Outrigger® II

Shoulder Immobilizer



The Outrigger II abducted shoulder immobilizer is now more comfortable and easier to apply.

Softer Neck/Shoulder Pad & Hand Strap

The new Outrigger II features a larger and softer neck/shoulder pad for improved comfort. The hand strap is also padded.

Adjustable Neck/Shoulder Pad

Unique fastening tab allows neck/shoulder pad to be easily repositioned on shoulder strap. Waist and shoulder straps also feature quick release buckles for easier application.

Breathable Sling

Sling envelope is now made entirely of air mesh spacer fabric for enhanced breathability and support.

Indications:

- Rotator cuff repairs
- Postoperative Bankhart procedures
- Capsular shifts
- Glenohumeral dislocation or subluxation
- Soft tissue shoulder strains or repairs



Pillow is designed to position the arm in approximately 15° of abduction.



Bilateral design fits left or right shoulder

Item #	Size	Measurements*
183752	Small	up to 13½"
183754	Medium	up to 14½"
183755	Large	up to 15½"
183756	X-Large	up to 17½"

*measure from Olecranon (tip of elbow) to MP joint (knuckles)

Includes exercise ball to promote circulation and provide hand rest.

Suggested HCPCS: L3670

