

# V-Strap™ Universal Wrist Support



When you want to provide your patients with the immobilization they need and the easy application and comfort they demand in order to be compliant, prescribe the V-Strap Wrist Support. The universal design fits a broad range of patient sizes and minimizes inventory costs.

### V-Strap Retention System

Provides ease of application since the "V" formed by the closure strap and the hook fastener prevents the strap from disengaging with the plastic buckle when applying or removing the wrist support.

### Superior Immobilization

The proximal and distal closure straps work in concert with the palmar and the dorsal aluminum stays to resist flexion and extension of the wrist. Palmar stay may be removed if desired.

### Universal Design

Dorsal stay casing may easily be repositioned to fit a wide range of patient sizes.

### Optimum Comfort

The polypropylene felt liner and outer nylon knit fabric allow moisture to evaporate away from the skin and enhance patient compliance. The anatomically contoured design and the adjustable thumb strap provide a comfortable fit and allow unrestricted motion of the fingers and thumb.

### Indications:

- Carpal tunnel syndrome
- Wrist sprains or strains
- Post cast removal
- Post operative care



8" wrist support



10 1/2" wrist & forearm support

black 8" left	black 8" right	black 10 1/2" left	black 10 1/2" right	size	wrist circ.
223450	223451	223460	223461	univ.	5 1/2" - 8 1/2"

Suggested HCPCS: L3908



rev 3/24

Medical Specialties Inc. 4600-K Lebanon Rd. Charlotte, NC 28227 p: 800-334-4143 f: 800-694-9060 p: 704-694-2434 f: 704-694-9060  
 e-mail: request@medspec.com www.medspec.com  
 © 2008-2024 Medical Specialties Inc.

