

Med Spec® SI Belt



With its unique design, the Med Spec SI Belt provides exceptional pelvic-sacral support. The superior materials of the Med Spec SI Belt allow for the brace to be worn comfortably underneath clothing while remaining in its proper position throughout daily activities. The Med Spec SI Belt provides your patients with the sacroiliac support they need and enhances their compliance with a comfortable and low profile design.

Features


- Outer inelastic nylon belt
- Posterior plastic panel & foam pad
- Two lateral SkinLoc™ pads
- Latex-free elastic liner

Benefits

- provides circumferential compression to help immobilize the sacroiliac joints
- provide additional circumferential compression to help immobilize the sacrum
- grip the skin or underlying clothing to limit the belt's migration
- allows for ease of application and is comfortable enough to be worn directly next to the skin



Item #	Size	Waist Dimen.
163302	small	26" - 31"
163304	medium	31" - 38"
163305	large	38" - 45"
163306	x-large	45" - 52"
163307	xx-large	52" - 59"



Posterior Plastic Panel & Foam Pad



May be worn directly next to the skin



Lateral SkinLoc Pad

PDAC Verified HCPCS: L0621

Indications:

- Sacroiliac joint dysfunction
- Sacroiliac joint pain



Medical Specialties Inc. 4600-K Lebanon Rd. Charlotte, NC 28227
 p: 800-334-4143 f: 800-694-9060 p: 704-694-2434 f: 704-694-9060
 e-mail: request@medspec.com www.medspec.com

© 2006-2024 Medical Specialties Inc.

