

# Supinator™

## PTT Stabilizer



The Supinator is a comfortable and relatively low-profile brace for the treatment of posterior tibial tendon dysfunction (PTTD). With its patent pending design, the Supinator reduces over-pronation by applying a supinating force.

### Arch Strap & Pad

The arch strap works in concert with a repositionable arch pad to lift the medial arch. The arch pad may also be removed if desired and to allow the Supinator to be worn with an orthotic.

### Heel Strap

The heel strap restricts eversion of the calcaneus (heel).

### Speed Lacer

Speed Lacer II closure system allows for ease of application and provides equal tension across all laces.

### Contraindications:

- Lateral ankle sprain
- Lateral ankle instability

### Indications:

- Posterior tibial tendon dysfunction (PTTD)
- Acquired adult flatfoot deformity
- Posterior tibial tendon insufficiency



left foot pictured above



right foot pictured above



left foot pictured above

Suggested HCPCS: L1902

Left	Right	Size	Ankle Circum.*
264251	264261	X-Small	10" - 11"
264252	264262	Small	11" - 12"
264254	264264	Medium	12" - 13"
264255	264265	Large	13" - 14"
264256	264266	X-Large	14" - 15"

\*measure ankle circumference according to diagram below



Men's Shoe Size	6-7	7 1/2-9	9 1/2-11	11 1/2-13	13 1/2-15
Narrow Foot		xs	s	m	l
Regular Foot	xs	s	m	l	xl
Wide Foot	s	m	l	xl	
Women's Shoe Size	5-6	6 1/2-8	8 1/2-10	10 1/2-12	12 1/2-14
Narrow Foot		xs	s	m	l
Regular Foot	xs	s	m	l	xl
Wide Foot	s	m	l	xl	

Medical Specialties Inc. 4600-K Lebanon Rd. Charlotte, NC 28227  
p: 800-334-4143 f: 800-694-9060 p: 704-694-2434 f: 704-694-9060  
e-mail: request@medspec.com www.medspec.com Patent Pending

Supinator is a trademark of Medical Specialties Inc. © 2018-2024 Medical Specialties Inc.



rev 9/21