

Ryno Lacer® II

Wrist & Thumb Support



With its patent pending design, the new Ryno Lacer II provides equal tension across all laces offering supreme comfort, ease of application and superior immobilization of the wrist and thumb.

Lacer II Closure System

Unique lace closure design provides equal tension across all laces to easily generate uniform compression for unparalleled support. The tension can be easily adjusted for patients with disproportionate forearm and wrist circumferences.

Comfort and Compliance

New softer proprietary felt liner contains polypropylene and wicks moisture away from the skin. Memory foam padding in palm and around base of thumb enhances comfort and support.

Radial, Palmar and Dorsal Stays


Three aluminum stays comfortably immobilize the wrist and thumb. In addition, each stay is removable if desired.

Indications:

- DeQuervain's Tendonitis
- Gamekeeper's Thumb
- Post Cast Healing
- Select Scaphoid Injuries



Long Ryno Lacer II



Short Ryno Lacer II

Suggested HCPCS: L3807 or L3809

Long Black Left	Long Black Right	Size
223711	223721	x-small
223712	223722	small
223714	223724	medium
223715	223725	large
223716	223726	x-large

Sizing Chart	
Size	Wrist Circ.
x-small	4 ³ / ₄ " - 5 ¹ / ₂ "
small	5 ¹ / ₂ " - 6 ¹ / ₄ "
medium	6 ¹ / ₄ " - 7"
large	7" - 7 ³ / ₄ "
x-large	7 ³ / ₄ " - 8 ¹ / ₂ "

Short Black Left	Short Black Right	Size
223731	223741	x-small
223732	223742	small
223734	223744	medium
223735	223745	large
223736	223746	x-large



Medical Specialties Inc. 4600-K Lebanon Rd. Charlotte, NC 28227
 p: 800-334-4143 f: 800-694-9060 p: 704-694-2434 f: 704-694-9060
 e-mail: request@medspec.com www.medspec.com US Patent No. 11,883,316
 Ryno Lacer is a registered trademark of Medical Specialties Inc. © 2017-2024 Medical Specialties, Inc.

