

AKS™ Knee Support

with Metal Hinges & Straps



The AKS Knee Support offers your patients the unique ability to adjust compression and fit from top to bottom.



Adjustable Knee Sleeve

The distal and proximal openings allow for easy application and ability to fit various patient shapes. In addition, compression can be adjusted to the desired level.

Adjustable Patellar Buttress

Soft horseshoe buttress is repositionable and works in concert with central sleeve region to provide circumferential compression around the knee and focal compression around the patella.

Medial and Lateral Metal Hinges & Straps

Metal hinges and outer closure straps provide medial and lateral stability. Internal Skinloc™ pad and metal hinges function synergistically to resist distal migration and vertical bunching. Condylar pads enhance comfort.



Indications:

- Arthritis of the knee
- Bursitis of the knee
- Chondromalacia patella
- Lateral or medial collateral ligament sprains
- Lateral or medial instability of the knee
- Patellar tendinitis

Bilateral design fits left or right knee.

Neoprene version provides greater compression and warmth.

CoolFlex version is breathable and lightweight for superior comfort.

Suggested HCPCS:
L1820 or L1821

Neoprene	Size	Knee Circ.*	CoolFlex
117111	X-Small	10" - 12"	117121
117112	Small	12" - 14"	117122
117114	Medium	14" - 16"	117124
117115	Large	16" - 18"	117125
117116	X-Large	18" - 20"	117126
117117	XX-Large	20" - 22"	117127
117118	3X-Large	over 22"	117128

*measure circumference of knee at mid patella



rev 10/24

Medical Specialties Inc.
4600-K Lebanon Rd.
Charlotte, NC 28227
p: 800-334-4143 f: 800-694-9060
p: 704-694-2434 f: 704-694-9060
e-mail: request@medspec.com
www.medspec.com
Patent pending

©2010-2024 Medical Specialties Inc.