

Dynatrack™

patella stabilizer with **Hinges**



The DynaTrack's unique buttress design and sturdy hinges work in synergy to normalize patella tracking, provide medial and lateral support, and control knee hyperextension. All of the features incorporated into the DynaTrack with Hinges function in concert to provide your patients with greater support and the comfort they demand to ensure compliance.

Features

- Aluminum hinges
- Internal patella buttress
- Anterior closure design
- Internal SkinLoc pad
- Open popliteal

Benefits

- Offer superior medial and lateral stability
- Applies direct, adjustable pressure to patella and remains in position due to SkinLoc™ lining material
- Allows for ease of application and an optimum fit
- Grips proximal end of calf to inhibit distal migration
- Minimizes material bunching during flexion to improve comfort

Indications:

- Patella subluxation
- Patella dislocation
- Chronic patellar maltracking
- Patellofemoral pain syndrome
- Post lateral release
- Post patellar realignment



Left	Right	Size	Knee Circ.
117041	117051	x-small	10" - 12"
117042	117052	small	12" - 14"
117044	117054	medium	14" - 16"
117045	117055	large	16" - 18"
117046	117056	x-large	18" - 20"
117047	117057	xx-large	20" - 22"

*measure circumference of knee at mid patella

Suggested HCPCS: L1820 or L1821

Medical Specialties Inc. 4600-K Lebanon Rd. Charlotte, NC 28227
 p: 800-334-4143 f: 800-694-9060 p: 704-694-2434 f: 704-694-9060
 e-mail: request@medspec.com www.medspec.com
 U.S. Patent # 7,004,919 © 2009-2024 Medical Specialties Inc.

