

# Gripper™ ROM Hinged Knee Brace



The *Gripper*™ Range-of-Motion Hinged Knee Brace provides your patients with the stability and versatility they need while resisting distal migration.



16" Gripper ROM

### Control & Adjustability

Easy access polycentric Range-of-Motion hinges allows for quick positioning of flexion and extension stops at 0°, 15°, 30°, 60°, and 90°. ROM hinges have an additional padded cover to allow for athletic participation. Inelastic nylon straps work in concert with ROM hinges to control flexion and extension and provide medial and lateral stability.

### Compliance & Migration

Internal SkinLoc™ pad works in concert with the outer closure strap to grip the proximal end of the calf muscle and inhibit the brace from migrating distally, thereby greatly enhancing patient compliance.

### Durability & Comfort

ROM hinge pockets are reinforced on the inside with ballistic nylon material at wear points to prevent hinges from wearing through the pocket. Unique condylar pad design enhances comfort. Wraparound anterior closure allows for ease of application and an optimum fit.

### Indications:

- Meniscal injuries
- MCL or LCL injuries or sprains
- Medial or lateral knee instabilities

PDAC Verified HCPCS:

L1832 or L1833



12" Gripper ROM



The Gripper ROM Hinge is available in two different base materials: Neoprene and CoolFlex™. The neoprene version provides warmth and superior compression. The CoolFlex version is breathable and hypoallergenic. All versions of the Gripper are latex free.

16" Neoprene	16" CoolFlex	12" Neoprene	12" CoolFlex	Size	6" Above	Mid-Patella	6" Below
117132	117142	117162	117172	Small	13" - 16"	12" - 14"	12" - 14"
117134	117144	117164	117174	Medium	16" - 18"	14" - 16"	13" - 15"
117135	117145	117165	117175	Large	18" - 21"	16" - 18"	15" - 17"
117136	117146	117166	117176	X-Large	21" - 24"	18" - 20"	17" - 19"
117137	117147	117167	117177	XX-Large	24" - 27"	20" - 22"	19" - 21"
117138	117148			3X-Large	27" - 31"	22" - 24"	21" - 23"

\*Use mid-patella measurement for primary sizing info



Medical Specialties Inc. 4600-K Lebanon Rd. Charlotte, NC 28227 p: 800-334-4143 f: 800-694-9060  
 p: 704-694-2434 f: 704-694-9060 e-mail: request@medspec.com www.medspec.com  
 ©2010-2024 Medical Specialties Inc.

