



# Gripper™ 3

## 12" Hinged Knee Brace



The Gripper™ 3 Hinged Knee Brace offers you the ability to provide your patients with a custom fitted device that delivers superior medial & lateral stability and resistance to distal migration.



### Enhanced Support & Custom Fit

Multiple condylar pads may be repositioned medially or laterally to provide enhanced support and a custom fit for patients with a valgus or varus knee. Easy-peel outer closure straps may also be trimmed to provide a custom fit.

### Compliance & Migration

Internal SkinLoc™ pad works in concert with the outer closure strap to grip the proximal end of the calf muscle and inhibit the brace from migrating distally, thereby enhancing patient compliance.

### Durability

Metal hinge pockets are reinforced on the inside with a woven nylon fabric at wear points to prevent hinge from wearing through the pocket.

### Ease of Application & Comfort

The wraparound anterior closure allows for ease of application and an optimum fit. Four-way stretch fabric in the popliteal fossa minimizes bunching during flexion to improve comfort. Anterior and sides of brace are made with CoolFlex™ material that is breathable and supportive.

### Indications:

- Meniscal injuries
- MCL or LCL injuries or sprains
- Medial or lateral knee instabilities

Suggested HCPCS:  
L1820 or L1821

Item #	Size	6" Above	Mid-Patella	6" Below
117231	X-Small	11" - 13"	10" - 12"	10" - 12"
117232	Small	13" - 16"	12" - 14"	12" - 14"
117234	Medium	16" - 18"	14" - 16"	13" - 15"
117235	Large	18" - 21"	16" - 18"	15" - 17"
117236	X-Large	21" - 24"	18" - 20"	17" - 19"
117237	XX-Large	24" - 27"	20" - 22"	19" - 21"
117238	3X-Large	27" - 31"	22" - 24"	21" - 23"
117239	4X-Large	31" - 35"	24" - 26"	23" - 26"

\*Use mid-patella measurement for primary sizing info

Medical Specialties Inc. 4600-K Lebanon Rd. Charlotte, NC 28227 p: 800-334-4143 f: 800-694-9060  
p: 704-694-2434 f: 704-694-9060 e-mail: request@medspec.com www.medspec.com

©2024 Medical Specialties Inc.

