

# EpiGel<sup>®</sup> elbow support



A comfortable, low profile design that aids in the treatment of tennis elbow (lateral epicondylitis) or golfer's elbow (medial epicondylitis). The EpiGel elbow support provides relief with two adjustable gel pads that compress the extensor and flexor muscles without restricting circulation of the forearm. The soft gel pads also dampen vibration. The universal design fits a broad size range of patients.

## EpiGel

Item #	Size	Elbow Circ.
223410	univ	8" - 14"

\*measure forearm just below elbow

### Two gel pads

Compress extensor muscles to relieve lateral epicondylitis and flexor muscles to relieve medial epicondylitis without creating a tourniquet around the forearm. Gel pads also dampen vibration.

### Unique strap retention system

Allows for easy application.

### Universal design

Fits right or left arm. One size fits most.

### Low profile design

Less likely to catch on any object.

# EpiFoam<sup>™</sup> elbow support

## EpiFoam

Item #	Size	Elbow Circ.
223420	univ	8" - 14"

\*measure forearm just below elbow



### Single foam pad

Compress extensor muscles to relieve lateral epicondylitis or flexor muscles to relieve medial epicondylitis.

### Universal design

Fits right or left arm. One size fits most.

### Unique strap retention system

Allows for easy application.

### Low profile design

Less likely to catch on any object.

Medical Specialties Inc. 4600-K Lebanon Rd. Charlotte, NC 28227  
 p: 800-334-4143 f: 800-694-9060 p: 704-694-2434 f: 704-694-9060  
 e-mail: request@medspec.com www.medspec.com  
 EpiGel, EpiFoam and Med Spec are trademarks of Medical Specialties Inc.  
 © 2010-2024 Medical Specialties Inc.



Suggested HCPCS:  
A4467

