

Gripper™

Hinged Knee Brace 16"



The Gripper™ Hinged Knee Brace provides your athletes and patients with the medial and lateral stability they need while resisting distal migration.



Compliance & Migration

Internal SkinLoc™ pad works in concert with the outer closure strap to grip the proximal end of the calf muscle and inhibit the brace from migrating distally, thereby enhancing patient compliance

Durability

Straps are made of high tension elastic that allows for numerous applications without fraying or losing closure strength. Metal hinge pockets are reinforced on the inside with ballistic nylon material at wear points to prevent hinge from wearing through the pocket.

Ease of Application & Comfort

The wraparound anterior closure allows for ease of application and an optimum fit. Open popliteal minimizes material bunching during flexion to improve comfort.

Indications:

- Meniscal injuries
- MCL or LCL injuries or sprains
- Medial or lateral knee instabilities



Suggested HCPCS: L1820 or L1821

The Gripper 16" is available in two different base materials: Neoprene and CoolFlex™. The neoprene version provides warmth and superior compression. The CoolFlex version is breathable and hypoallergenic.



geared polycentric hinge

Neoprene Geared	CoolFlex Geared	Size	6" Above	Mid-Patella	6" Below
117072	117082	Small	13" - 16"	12" - 14"	12" - 14"
117074	117084	Medium	16" - 18"	14" - 16"	13" - 15"
117075	117085	Large	18" - 21"	16" - 18"	15" - 17"
117076	117086	X-Large	21" - 24"	18" - 20"	17" - 19"
117077	117087	XX-Large	24" - 27"	20" - 22"	19" - 21"

*Use mid-patella measurement for primary sizing info



Medical Specialties Inc. 4600-K Lebanon Rd. Charlotte, NC 28227 p: 800-334-4143 f: 800-694-9060
 p: 704-694-2434 f: 704-694-9060 e-mail: request@medspec.com www.medspec.com
 ©2010-2025 Medical Specialties Inc.

