



# Gripper™ 3

## Hinged Knee Brace 12"



The Gripper™ 3 Hinged Knee Brace offers you the ability to provide your patients with a custom fitted device that delivers superior medial & lateral stability and resistance to distal migration.



### Enhanced Support & Custom Fit

Multiple condylar pads may be repositioned medially or laterally to provide enhanced support and a custom fit for patients with a valgus or varus knee. Easy-peel outer closure straps may also be trimmed to provide a custom fit.

### Compliance & Migration

Internal SkinLoc™ pad works in concert with the outer closure strap to grip the proximal end of the calf muscle and inhibit the brace from migrating distally, thereby enhancing patient compliance.

### Durability

Metal hinge pockets are reinforced on the inside with a woven nylon fabric at wear points to prevent hinge from wearing through the pocket.

### Ease of Application & Comfort

The wraparound anterior closure allows for ease of application and an optimum fit. Four-way stretch fabric in the popliteal fossa minimizes bunching during flexion to improve comfort. Anterior and sides of brace are made with CoolFlex™ material that is breathable and supportive.

### Indications:

- Meniscal injuries
- MCL or LCL injuries or sprains
- Medial or lateral knee instabilities

Suggested HCPCS:  
L1820 or L1821

| Item # | Size     | 6" Above  | Mid-Patella | 6" Below  |
|--------|----------|-----------|-------------|-----------|
| 117231 | X-Small  | 11" - 13" | 10" - 12"   | 10" - 12" |
| 117232 | Small    | 13" - 16" | 12" - 14"   | 12" - 14" |
| 117234 | Medium   | 16" - 18" | 14" - 16"   | 13" - 15" |
| 117235 | Large    | 18" - 21" | 16" - 18"   | 15" - 17" |
| 117236 | X-Large  | 21" - 24" | 18" - 20"   | 17" - 19" |
| 117237 | XX-Large | 24" - 27" | 20" - 22"   | 19" - 21" |
| 117238 | 3X-Large | 27" - 31" | 22" - 24"   | 21" - 23" |
| 117239 | 4X-Large | 31" - 35" | 24" - 26"   | 23" - 26" |

\*Use mid-patella measurement for primary sizing info

Medical Specialties Inc. 4600-K Lebanon Rd. Charlotte, NC 28227 p: 800-334-4143 f: 800-694-9060  
p: 704-694-2434 f: 704-694-9060 e-mail: request@medspec.com www.medspec.com

©2025 Medical Specialties Inc.

