

Supinator™

PTT Stabilizer



The Supinator is a comfortable, durable and relatively low-profile brace for the treatment of posterior tibial tendon dysfunction (PTTD) or progressive collapsing foot deformity (PCFD).

Midfoot Strap

The midfoot strap works in concert with a repositionable arch pad to lift the medial arch. The arch pad may also be removed if desired or to allow the Supinator to be worn with an orthotic.

Hindfoot Strap

The hindfoot strap restricts eversion of the calcaneus (heel).

Speed Lacer

Speed Lacer II closure system allows for ease of application and provides equal tension across all laces.

Contraindications:

- Lateral ankle sprain
- Lateral ankle instability

Indications:

- Posterior tibial tendon dysfunction (PTTD)
- Progressive collapsing foot deformity (PCFD)
- Acquired adult flatfoot deformity (AAFD)



right foot pictured above



left foot pictured above

Suggested HCPCS: L1902

Left	Right	Size	Ankle Circum.*
264251	264261	X-Small	10" - 11"
264252	264262	Small	11" - 12"
264254	264264	Medium	12" - 13"
264255	264265	Large	13" - 14"
264256	264266	X-Large	14" - 15"

*measure ankle circumference according to diagram below



left foot pictured above

Men's Shoe Size	6-7	7 1/2-9	9 1/2-11	11 1/2-13	13 1/2-15
Narrow Foot		xs	s	m	l
Regular Foot	xs	s	m	l	xl
Wide Foot	s	m	l	xl	
Women's Shoe Size	5-6	6 1/2-8	8 1/2-10	10 1/2-12	12 1/2-14
Narrow Foot		xs	s	m	l
Regular Foot	xs	s	m	l	xl
Wide Foot	s	m	l	xl	

Medical Specialties Inc. 4600-K Lebanon Rd. Charlotte, NC 28227
 p: 800-334-4143 f: 800-694-9060 p: 704-694-2434 f: 704-694-9060
 e-mail: request@medspec.com www.medspec.com Patent Pending

Supinator is a trademark of Medical Specialties Inc. © 2018-2025 Medical Specialties Inc.



rev 9/25