

# Suede Thumb Support



Want to provide your patients with an easy to apply, comfortable, and low profile splint for deQuervain's Tendonitis? Then, the Med Spec **Suede Thumb Support** may be the solution you are looking for.

## Indications:

- deQuervain's Tendonitis
- Gamekeeper's Thumb
- Thumb sprains or strains
- Post cast removal
- Post operative care

## Optimum Comfort & Support

The perforated suede and polypropylene felt liner allow for moisture to evaporate away from the skin. In addition, the anatomically contoured aluminum stay and soft elastic closure straps comfortably immobilize the MP and CMC joints of the thumb.

## V-Strap Retention System

Provides extraordinary ease of application since the "V" formed by the elastic closure strap and the hook fastener prevents the strap from disengaging with the plastic buckle when applying or removing the thumb support.

## Dual Thumb Closure Straps


The inner elastic strap and the outer closure strap work in concert to allow for an easy and secure closure around the thumb.



Short Suede Thumb Support



Long Suede Thumb Support



Short	Size	Long
223180	universal	223190

PDAC Verified HCPCS: L3923

Medical Specialties Inc.  
4600-K Lebanon Rd.  
Charlotte, NC 28227

p: 800-334-4143 f: 800-694-9060  
p: 704-694-2434 f: 704-694-9060  
e-mail: request@medspec.com  
www.medspec.com

© 2007-2026 Medical Specialties Inc.



rev 1/26