

Med Spec[®]

Boxer Splint II



Treat Multiple Fractures With One Splint

Optimize fracture management with a single, versatile splint — ideal for practices with limited DME inventory space.



Radial Gutter

Versatile Design

Functions as a Boxer's Fracture Splint, Radial Gutter Splint, and 3rd & 4th Metacarpal Buddy Splint.

Features malleable aluminum stays that are easily contoured for a customized fit and precise immobilization of the MP joints in the desired degree of flexion.

Comfort & Compliance

Soft-lined shell fabric provides breathable padding to enhance all-day comfort.

Memory-foam reinforcement enhances comfort by reducing pressure points and irritation.

Ease of Application

Integrated inner finger sleeve simplifies application and helps maintain alignment.

Sleeve minimizes skin-on-skin contact to help reduce maceration risk.

Finger sleeve may be trimmed if desired to achieve a custom fit.



Boxer's Fracture



3rd & 4th Buddy Splint

Indications:

- Metacarpal fractures
- Phalangeal injuries
- MP (MCIP) joint injuries
- PIP joint injuries

Suggested HCPCS:
L3807 or L3809

Part ID	Size	Palm Circ.*
224571	XS LT-Ulnar/RT-Radial	under 6"
224574	SM/MD LT-Ulnar/RT-Radial	6" - 9"
224576	LG/XL LT-Ulnar/RT-Radial	over 9"
224581	XS RT-Ulnar/LT-Radial	under 6"
224584	SM/MD RT-Ulnar/LT-Radial	6" - 9"
224586	LG/XL RT-Ulnar/LT-Radial	over 9"

*Measure circumference at palmar crease



Medical Specialties Inc. 4600-K Lebanon Rd. Charlotte, NC 28227 p: 800-334-4143 f: 800-694-9060
p: 704-694-2434 f: 704-694-9060 e-mail: request@medspec.com www.medspec.com
©2025-2026 Medical Specialties Inc.

